

NEW: Fractionation on DEXTech Instruments Saves Time and Money



One of the many strong points of our DEXTech series of instruments for the clean-up of PCDD/Fs and PCBs is the flexibility in fractionation possibilities. Just by pushing a button, customers can choose between different methods with different fraction options. Now, on multiple customer request, we have developed a fractionation option on the DEXTech instruments where all dioxins like-PCBs (DL-PCBs) are fractionated with the PCDD/Fs.

Sometimes labs are only interested in the results of the DL-PCBs and the PCDD/Fs. Conventional clean-up systems can only fractionate the non-ortho PCBs with the PCDD/Fs separated from the other mono-ortho PCBs. This means the analytical lab has to do two analytical runs on the GC to determine the result for DL-PCBs and PCDD/Fs. With the new DEXTech clean-up fractionation method, the DL-PCBs now can be measured together with PCDD/Fs in one GC run. As GC-runtime is one of the major bottlenecks in the analysis of PCDD/Fs and PCBs having only one GC run for the determination of PCDD/Fs and DL-PCBs saves a lot of time and money. For detailed information about the recoveries, please see next page.

Alox PLUS Method

Fraction 1 (24 mL*):
Mono-ortho-PCB
+ ndI-PCB + PBDE

Fraction 2 (10 mL):**
Non-ortho-PCB
+ PCDD/F
+ PCNs

Total process time: 65 min
Total process time SMART: 45 min

Also available as DCM Free Method!

Alox PURE Method

Fraction 1 (24 mL*):
Mono-ortho-PCB
+ ndI-PCB
+ non-ortho-PCB

Fraction 2 (10 mL):**
PCDD/F

Total process time: 72 min
Total process time SMART: 52 min

*UNIQUE FRACTIONATION:
matching best the analytical
setup in many laboratories*

NEW

DL - PCB + PCDD/F Method

Fraction 1 (24 mL*):
+ ndI-PCB

Fraction 2 (10 mL):**
DL-PCB+PCDD/F

Total process time: 65 min
Total process time SMART: 45 min

* Dichloromethane / n-hexane ** Toluene

For further information, please contact us using the details below.



		1018-4 Standard carbon column 4g Olive oil	1018-5 New carbon column 4g Olive oil	1018-6 New carbon column 4g Olive oil	1018-2 New carbon column 4g Fish oil	1018-3 New carbon column 4g Fish oil					
		F2	F2	F2	F2	F2					
13C- recoveries [%]	13C-2,3,7,8-TCDF	72	75	73	73	76					
	13C-1,2,3,7,8-PeCDF	72	76	74	83	83					
	13C-2,3,4,7,8-PeCDF	70	75	73	89	88					
	13C-1,2,3,4,7,8-HxCDF	87	77	90	88	85					
	13C-1,2,3,6,7,8-HxCDF	91	77	91	94	91					
	13C-2,3,4,6,7,8-HxCDF	90	79	89	93	88					
	13C-1,2,3,7,8,9-HxCDF	94	85	94	98	94					
	13C-1,2,3,4,6,7,8-HpCDF	98	87	97	98	94					
	13C-1,2,3,4,7,8,9-HpCDF	92	82	90	96	92					
	13C-1,2,3,4,6,7,8,9-OCDF	94	86	88	92	85					
	13C-2,3,7,8-TCDD	70	73	71	82	75					
	13C-1,2,3,7,8-PeCDD	78	76	71	82	89					
	13C-1,2,3,4,7,8 -HxCDD	104	94	109	105	103					
	13C-1,2,3,6,7,8-HxCDD	79	74	82	85	79					
	13C-1,2,3,7,8,9-HxCDD	101	88	101	99	94					
	13C-1,2,3,4,6,7,8-HpCDD	90	87	91	95	90					
	13C-1,2,3,4,6,7,8,9-OCDD	83	77	80	86	79					
	PCB-#81L	81	77	71	81	76					
	PCB-#77L	81	75	70	82	80					
	PCB-#126L	107	90	83	115	114					
PCB-#169L	109	93	81	127	121						
		F1	F2	F1	F2	F1	F2	F1	F2	F1	F2
13C- recoveries [%]	PCB-#28L	86	3	26	62	20	64	12	80	31	62
	PCB-#52L	70	1	95	1	91	1	87	1	92	1
	PCB-#101L	88	1	85	2	84	2	79	1	83	1
	PCB-#153L	92	1	87	2	85	2	85	1	82	1
	PCB-#138L	97	1	92	3	89	3	95	3	99	3
	PCB-#180L	94	1	83	4	81	3	82	2	84	2
	PCB-#81L	0	91	0	90	0	89	0	95	0	97
	PCB-#77L	0	86	0	86	0	84	0	93	0	94
	PCB-#126L	0	98	0	96	0	93	0	99	0	99
	PCB-#169L	1	98	2	101	2	94	0	92	2	91
	PCB-#123L	80	8	0	84	0	82	0	82	0	81
	PCB-#118L	91	10	0	82	0	81	0	88	0	86
	PCB-#114L	86	4	3	79	2	77	1	81	4	76
	PCB-#105L	80	9	0	80	0	76	0	84	0	83
	PCB-#167L	79	14	0	87	0	84	0	86	0	85
	PCB-#156L	83	10	0	85	0	80	0	83	0	83
	PCB-#157L	80	12	2	86	2	81	1	83	1	83
	PCB-#189L	84	12	1	84	1	79	1	87	1	84